



# The hemostat

Easy. Effective.

mediCipio<sup>®</sup>A

# mediCipio<sup>A</sup>

The patented two-layer structure of mediCipio A combines the excellent hemostyptic effect, absorbency and adhesiveness of oxidized, regenerated cellulose (ORC) with the stability and compatibility of porcine collagen, without the use of human fibrinogen and thrombin.

**Easy.**

To handle.


**Effective.**

To achieve hemostasis.




## Easy handling due to best material

- Ready to use - without prior moistening.
- Easy and flexible adaptation of the soft fleece to the wound (Fig. 1).
- mediCipio A sticks to the wound, not to the glove or instruments.
- Uncomplicated due to the absence of human thrombin: no coagulation in vessels or during autotransfusion.



Malleable already as a dry product



Excellent visibility for postoperative inspection


- On contact with tissue fluid, the matrix does not discolour, but becomes semi-transparent and allows the assessment of hemostasis (Fig. 2).
- Available in different sizes and easily cut to size sterilely if required.
- Suitable for use in the trocar.



## Ensuring medical success effectively

- Safe and reliable hemostasis
- Promotion of wound healing and best compatibility through porcine collagen
- No increase in volume – therefore also applicable near nerves
- Long term stability – no detachment of fragments of the fleece
- No allergenic additives, no risks from human thrombin or fibrinogen
- mediCipio A is completely absorbed within a few weeks

*Surgeons accomplish amazing medical achievements – the hemostat should keep up.*



Evaluation of the amount of fluid in the redon-drainage system after thyroidectomy

## Proven effectiveness: Best hemostyptic quality

In a prospective and randomized clinical trial with 150 patients it was confirmed that mediCipio A is as effective in inducing hemostasis as a product based on human thrombin and fibrinogen. Safe hemostasis was also confirmed when platelet aggregation inhibitors were used.

## Easy medical device – simply advantageous

- Medical device – restrictions for medicinal products do not apply
- Attractive pricing
- Cost savings through process optimization
- Storable at room temperature
- Made in Germany with strict quality control from raw materials to the final product
- Available in various sizes

Size	Content	PPN
1 x 5 cm	5 pcs.	11 1417 4850 60
2,5 x 3 cm	5 pcs.	11 1417 4867 50
5 x 5 cm	2 pcs.	11 1417 4873 16
5 x 10 cm	2 pcs.	11 1417 4904 63

Each individually sterile packed

Can be cut to size under sterile conditions

**Further information and orders:**

**Tel.: +49(0) 37190 96 83 8**



**medichema**

**mediCipio**

is a brand of

**medichema GmbH**

Weststraße 57

09112 Chemnitz

Tel.: +49(0) 371 90 96 83 8

Fax: +49(0) 371 90 96 83 9

[info@medichema.com](mailto:info@medichema.com)

[www.medichema.com](http://www.medichema.com)